



|   |   |         |
|---|---|---------|
|  | <b>ANALYSENZERTIFIKAT<br/>CERTIFICATE OF ANALYSIS</b> | DQ6MBA0 |
|---|---|---------|

**TECHNOCLOT® PT Owren manual 10x10mL**

REF      5005037  
LOT      6M72BB0.01  
      30-Apr-2019

| Analytische Prozedur<br>Analytical Procedure                         |   | Anforderung<br>Requirement        | Ergebnis<br>Result |    |
|--|---|-----------------------------------|--------------------|----|
| Verwendbarkeit<br>Usability<br>KC 10                                 | CCN<br><span style="border: 1px solid black; padding: 1px;">LOT</span> 1P37B01  | 70.3 – 201.8 %<br>0.82 – 1.16 INR | 81.5<br>1.06       | OK |
|  | CCAK<br><span style="border: 1px solid black; padding: 1px;">LOT</span> 2B35B00 | 13.1 – 19.8 %<br>2.64 – 3.71 INR  | 18.6<br>2.82       | OK |
| Etiketten und Barcode Überprüfung<br>Label and bar code Verification |   |                                   | OK                 |    |

**Durchführung von  
Prepared by**

*Lucas*  
**Lucas ROSENMAYR**

**10. MAI 2017**

Datum / Date:

**Freigabe von  
Approved by**

*Stefan*  
**Stefan STEINER**

**10. MAI 2017**

Datum / Date: