

# Thrombodynamics Analyser

## Thrombodynamics assay



26.04.2018 09:01

Measurement title: Normal Coagulation TD  
 Plasma type: Fresh PFP  
 Sample code: Example TDX #1

**TDX**



Fibrin dynamics parameters	Unit	Parameter	Value	Norm.
Rate of clot growth	um/min	V	26,0	20,0-29,0
Lag time	min	Tlag	1,0	0,6-1,5
Initial rate of clot growth	um/min	Vi	50,2	38,0-56,0
Stationary rate of clot growth	um/min	Vst	26,0	20,0-29,0
Clot size at 30 min	um	CS	1122	800-1200
Clot density	a.u.	D	28271	15000-32000
Spontaneous clotting	min	Tsp		30,0-



Comment:

Tech info:

**ATTENTION: Research use only. Measurement results shouldn't be used for diagnostic purposes.**

S/N: 000000, Operator: Operator, report created 29.03.2019 12:12