

# APC™ RESISTANCE TESTING

Gold Standard for FV Leiden Screening

See reverse for details & ordering information.

**The Diapharma team's deep knowledge base translates into over 500 solutions for your hemostasis needs.**

**Questions about any product in our lineup?  
These professionals are just a click or phone call away.**



diapharma

800.526.5224

info@diapharma.com

# APC™ RESISTANCE TESTING

A time- and money-saving alternative to individual genetic testing, Coatest® APC™ Resistance testing is an ideal solution for physicians and labs. Each kit was designed to be easy to use for in-house screening, offering fast results and low costs for both large- and small-volume labs. The most common FV Leiden screening test performed, it offers unmatched sensitivity for the FV:Q506 mutation – close to 100% specificity – and is applicable to patients on heparin or warfarin.



## APC™ RESISTANCE TESTING KIT FEATURES

IDEAL FOR:

- Labs looking to bring testing in-house
- Labs looking for cost savings
- Ordering physicians looking for a sensitive assay with rapid turnaround
- Saves the lab from having to perform genetic test on all patients
- Economical alternative to genetic testing

## DESCRIPTION

Coatest® APC™ Resistance V  
FV Leiden screening for larger volume labs

Coatest® APC™ Resistance VS  
FV Leiden screening for smaller volume labs

Coatest® APC™ Resistance  
For APC Resistance phenotype testing

Coatest® APC™ Resistance  
Control Plasma Level 1  
Normal Control Plasma

Coatest® APC™ Resistance  
Control Plasma Level 2  
Heterozygous Control Plasma

## REFERENCE & PACKAGING

K823120

K823138

K822643

C822650 • 5 x 1mL

C822668 • 5 x 1mL

800.526.5224      info@diapharma.com

ML-00-00118Rev04

