

DiaPharma Streptokinase

Data Sheet

Catalog #: B811111

Item Description: Streptokinase (Approximately 0.88 MIU)

Lot #: B-121901

Re-assay Date: 12/31/2021

All units of plasma have been tested for absence of HB_sAg (third generation test), anti-HIV (FDA licensed test) and found to be negative. This material should be handled as if capable of transmitting disease.

CAUTION:

NOT SUITABLE OR INTENDED FOR USE IN HUMANS OR FOOD ANIMALS

NOT FOR STERILE USE OR INJECTION.

IN VITRO USE ONLY.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product: Streptokinase (purified preparation of a bacterial protein elaborated by group C b-hemolytic streptococci). Supplied as white, water-soluble, lyophilized powder (non-sterile)

Material contains 11 mg sodium L-glutamate and 13.0 mg albumin (human) per 1,000,000 IU (1 MIU) of Streptokinase.

Contains no preservatives.

Amount:

per vial: Approximately **880,090 IU** (0.88 MIU) per vial

Prepared By: Jennifer J. Kiblinger **Date:** 12/20/2019

Printed Name: Jennifer J. Kiblinger

Title: Scientific Director