

COAMATIC[®] HEPARIN

Anti-Xa Testing

See reverse for details & ordering information.

The Diapharma team's deep knowledge base translates into over 500 solutions for your hemostasis needs.

**Questions about any product in our lineup?
These professionals are just a click or phone call away.**



diapharma

800.526.5224

info@diapharma.com

COAMATIC® HEPARIN

Coamatic® Heparin measures the ability of heparin to catalyze the inhibition of FXa by antithrombin. The anti-Xa assay is the method of choice when monitoring LMW heparin therapy and is also suitable as a replacement for the APTT test when monitoring therapy with UF Heparin.



ASSAY FEATURES

- Simplify Heparin testing
- Undiluted assay with superior reagent stability

DESCRIPTION

Chromogenix® Coamatic® Heparin

ALSO AVAILABLE

Chromogenix® Coatest® Heparin

Chromogenix® LMW Heparin Calibration Plasma

Chromogenix® LMW Heparin Control Plasma

REFERENCE & PACKAGING

K823393 • Chromogenic Kit/200 tests

REFERENCE & PACKAGING

K255539 • Chromogenic Kit/400 tests

C823500 • 4x3x1ml

C823492 • 4x2x1ml

800.526.5224 info@diapharma.com

ML-00-00183Rev01

