



## BRINGING CHROMOGENIC FACTOR IX DETERMINATION INTO THE NEW ERA

Rossix, a specialized company in hemostasis product development, presents **Rox Factor IX (FIX)** and **Rox FIX-A (FIXa)**, exclusively available from Diapharma in the US and Canada.




# MEASURING FACTOR IX AND FACTOR IX<sub>a</sub> HAS NEVER BEEN EASIER

## Rox Factor IX

- Chromogenic assay for Factor IX activity in plasma or concentrates
- High sensitivity at low levels, as little as 1% FIX in plasma
- Study of new FIX variants
  - Accurate pharmacokinetics in hemophilia B plasma samples
  - Does not overestimate FIX activity compared to certain aPTT reagents in quantifying FIX variants
- Reliable potency assignments of recombinant and plasma-derived FIX
- No interference from LA
- FIX deficient plasma not required
- Assay time less than 20 minutes
- FIX activity in animal models

## Rox Factor IX<sub>a</sub>

- Determine pre-activated FIX in FIX concentrates
- Pharmacopeial assay
- FIX deficient plasma not required
- Highly sensitive and robust. Measure linearly down to 0.02 mIU FIX<sub>a</sub> per U of FIX



### Ordering Information

Rox Factor IX	Chromogenic	100 determinations	900020
Rox FIX-A	Chromogenic	100 determinations	950030

Customer Service 800.526.5224  
 Technical Support 800.447.3846  
 diapharma.com

ML-00-00009Rev02

Rosix products are for research use only in the US and Canada

