

TECHNOZYM® ADAMTS-13
Fluorogenic ELISA

See reverse for details & ordering information.

The Diapharma team's deep knowledge base translates
into over 500 solutions for your hemostasis needs.

Questions about any product in our lineup?
These professionals are just a click or phone call away.



diapharma


800.526.5224

info@diapharma.com

TECHNOZYM® ADAMTS-13

Fluorogenic ELISA

TECHNOZYM® ADAMTS-13 Activity/Antigen Kit is a fluorogenic ELISA for the determination of ADAMTS-13 activity and antigen in plasma in research methods. ADAMTS-13 is the enzyme that cleaves vWF under shear stress condition. Dysfunction of this enzyme leads to unusually large vWF multimers and increased platelet aggregation.



ASSAY FEATURES

- Allows determination of ADAMTS-13 activity below 5%
- Samples used undiluted
- Excel evaluation software available

DESCRIPTION

12 / 6 x 8 Anti-ADAMTS-13 Strips
1 x 0.3 ml Anti-ADAMTS-13 Conjugate
1 x 80 ml Wash Buffer Concentrate
1 x 90 ml Incubation Buffer
2 / 1 x 3 ml Activity Substrate
1 x 6 ml Antigen Substrate
1 x 0.7 ml Stable Peroxidase Solution
1 x 6 ml Stop Solution
7 x 0.5 ml Calibrators and Controls

REFERENCE & PACKAGING

5450501 • 2 x 48 Tests
5450551 • 48 Tests

ADDITIONAL CALIBRATORS & CONTROLS

Additional Calibrator Set for TECHNOZYM® ADAMTS-13 Activity and Antigen
5 x 0.5 ml Calibrators

REFERENCE & PACKAGING

5450561

Additional High and Low Control for TECHNOZYM® ADAMTS-13 Activity and Antigen

5450563

800.526.5224 info@diapharma.com

RESEARCH USE ONLY • ML-00-00152REV01

