

TECHNOZYM® ADAMTS-13
Inhibitor ELISA

See reverse for details & ordering information.

The Diapharma team's deep knowledge base translates
into over 500 solutions for your hemostasis needs.

Questions about any product in our lineup?
These professionals are just a click or phone call away.




diapharma

800.526.5224

info@diapharma.com

TECHNOZYM® ADAMTS-13 Inhibitor ELISA

TECHNOZYM® ADAMTS-13 INH is an ELISA for the determination of antibodies directed against ADAMTS-13 in research methods.



ASSAY FEATURES

- Ready-to-use coated ELISA plate
- Low background
- No false positive results with samples containing high amounts of IgG in general
- Specific for anti-ADAMTS-13 IgG
- Reproducible and simple – easy to standardize
- Fast and easy to perform

DESCRIPTION

12 x 8 Recombinant ADAMTS-13 Strips
1 x 0.3 ml Anti-IgG Conjugate
1 x 80 ml Wash Buffer Concentrate
1 x 90 ml Incubation Buffer
1 x 12 ml TMB Substrate
1 x 12 ml Stop Solution
7 x 0.5 ml Calibrators and Controls

REFERENCE & PACKAGING

5450401 • 96 Tests
5450451 • 48 Tests

ADDITIONAL CALIBRATORS & CONTROLS

Additional Calibrator Set for
TECHNOZYM® ADAMTS-13 INH

5450461 • 5 x 0.5 ml

Additional High and Low Control for
TECHNOZYM® ADAMTS-13 INH

5450463 • 2 x 0.5 ml

800.526.5224 info@diapharma.com

RESEARCH USE ONLY • ML-00-00153REV01

