

TECHNOZYM® ADAMTS-13

Antigen ELISA

See reverse for details & ordering information.

The Diapharma team's deep knowledge base translates into over 500 solutions for your hemostasis needs.

Questions about any product in our lineup?
These professionals are just a click or phone call away.



diapharma


800.526.5224

info@diapharma.com

TECHNOZYM® ADAMTS-13

Antigen ELISA

TECHNOZYM® ADAMTS-13 Antigen assay is an ELISA for the determination of ADAMTS-13 Antigen concentration in plasma in research methods. ADAMTS-13 is an enzyme that specifically cleaves vWF under high shear stress.



ASSAY FEATURES

- Detection of low ADAMTS-13 Antigen concentrations
- No need for any special technique or instrumentation, just **standard ELISA equipment**

DESCRIPTION

12 x 8 Anti-ADAMTS-13 Strips
1 x 0.3 ml Anti-ADAMTS-13 Conjugate
1 x 80 ml Wash Buffer Concentrate
1 x 90 ml Incubation Buffer
1 x 12 ml TMB Substrate
1 x 12 ml Stop Solution
7 x 0.5 ml Calibrators and Controls

REFERENCE & PACKAGING

5450601 • 96 Tests

ADDITIONAL CALIBRATORS & CONTROLS

Additional Calibrator Set for
TECHNOZYM® ADAMTS-13 Antigen

REFERENCE & PACKAGING

5450661 • 5 x 0.5 ml

Additional High and Low Control for
TECHNOZYM® ADAMTS-13 Antigen

5450663 • 2 x 0.5 ml

800.526.5224 info@diapharma.com

RESEARCH USE ONLY • ML-00-00154REV01

