

TECHNOCHROM® C1-INH

See reverse for details & ordering information.

The Diapharma team's deep knowledge base translates into over 500 solutions for your hemostasis needs.

Questions about any product in our lineup?
These professionals are just a click or phone call away.




diapharma

800.526.5224

info@diapharma.com

TECHNOCHROM® C1-INH

TECHNOCHROM® C1-Inhibitor is a reagent kit for the chromogenic determination of C1-esterase inhibitor (C1-INH) in research methods. C1-INH deficiency leads to hereditary angioedema.



TECHNOCHROM® C1-INH is a simple and fast assay for determining C1 activity in plasma or concentrates.

DESCRIPTION

- 1 x 3 ml Substrate
- 1 x 3 ml C1-Esterase
- 1 x 25 ml Sample Buffer A
- 1 x 20 ml Reaction Buffer B
- 1 x 1 ml Reference Plasma
- 2 x 1 ml Normal and Abnormal Controls

REFERENCE

5345003

PACKAGE

~ 30 Tests

To order TECHNOCHROM® C1-INH or for any questions about our entire product line, give us a call or visit us online.

800.526.5224 info@diapharma.com

RESEARCH USE ONLY • ML-00-00155REV01

