

TECHNOZYM® ADAMTS-13  
Activity ELISA

See reverse for details & ordering information.

The Diapharma team's deep knowledge base translates  
into over 500 solutions for your hemostasis needs.

Questions about any product in our lineup?  
These professionals are just a click or phone call away.



diapharma


800.526.5224

info@diapharma.com

# TECHNOZYM® ADAMTS-13

Activity ELISA

TECHNOZYM® ADAMTS-13 Activity assay is an ELISA for the determination of ADAMTS-13 activity in plasma in research methods. ADAMTS-13 is an enzyme that specifically cleaves unusually large von Willebrand Factor multimers, which include platelet thrombus formation under high shear stress.



## ASSAY FEATURES

- Severe ADAMTS-13 deficiency detection – sensitivity down to 0.5%
- No need for any special technique or instrumentation, just **standard ELISA equipment**
- No interferences from bilirubin or hemoglobin

## DESCRIPTION

12 x 8 Anti-GST ELISA Test Strips  
1 x 12 ml HRP Conjugate  
1 x 30 ml Reaction Buffer  
1 x 53 ml Wash Buffer Concentrate  
1 x 12 ml TMB Substrate  
2 x 6 ml GST-vWF73 Substrate  
1 x 12 ml Stop Solution  
8 x 0.5 ml Calibrators and Controls

## REFERENCE & PACKAGING

5450701 • 96 Tests

## ADDITIONAL CALIBRATORS & CONTROLS

Additional Calibrator Set for  
TECHNOZYM® ADAMTS-13 Activity

## REFERENCE & PACKAGING

5450761 • 6 x 0.5 ml

Additional High and Low Control for  
TECHNOZYM® ADAMTS-13 Activity

5450763 • 2 x 0.5 ml

800.526.5224 info@diapharma.com

RESEARCH USE ONLY • ML-00-00151REV01

