

C	PRODUCT NAME	UNITS	REFERENCE VALUE			WHO STD
			ACL Classic	ACL Elite/Elite Pro/8/9/10000	ACL TOP Family / ACL TOP Family 50 Series†	
Fibrinogen	PT-Fibrinogen (1)	mg/dL g/L				
	PT-Fibrinogen HS PLUS (1)	mg/dL g/L				
	PT-Fibrinogen Recombinant (1)	mg/dL g/L			NA	
	RecombiPlasTin 2G (1)	mg/dL g/L				
	ReadiPlasTin (1)	mg/dL g/L	NA			
	Fibrinogen-C	mg/dL g/L				
	Q.F.A.	mg/dL g/L	NA			
Factor II	Factor II Deficient Plasma	%				
Factor V	Factor V Deficient Plasma	%				
Factor X	Factor X Deficient Plasma	%				
Factor VII	Factor VII Deficient Plasma	%				
Factor VIII	Factor VIII Deficient Plasma	%				
Factor VIII	Chromogenic Factor VIII	%	NA		NA	
Factor IX	Factor IX Deficient Plasma	%				
Factor XI	Factor XI Deficient Plasma	%				
Factor XII	Factor XII Deficient Plasma	%				
AT	Antithrombin	%				
PLG	Plasminogen	%				
PI	Plasmin Inhibitor	%				
PC	ProClot	%				
	Protein C	%				
PS	ProS (Clotting-functional)	%				
	Protein S Activity* (Clotting-functional)	%				
	Free Protein S (Latex-immunological)	%	NA			
VWF:Ag	von Willebrand factor Antigen	%		**		
VWF Activity	von Willebrand factor Activity	%	NA			
VWF:RCo	Von Willebrand Factor Ristocetin Cofactor Activity*	%	NA	NA		
Factor XIII	Factor XIII Antigen	%	NA	NA		

Note / Nota / Note

*Not available in all countries/No disponible en todos los países/Non disponible dans certains pays

**ACL 6000/7000

(1) PT based
Basado en TP
Derivé du TP

† ACL TOP Family 50 Series = ACL TOP 350 CTS; ACL TOP 550 CTS; ACL TOP 750; ACL TOP 750 CTS; ACL TOP 750 LAS

LOT _____



Symbols used / Símbolos utilizados / Symboles utilisés

UNITS	Measurement units / Unidades de medida / Unités de mesure
REFERENCE VALUE	Reference value / Valor de referencia / Valeurs de référence
WHO STD	WHO Standard Code / Código del Estándar OMS / Code standard OMS
PRODUCT NAME	IL Product name / Nombre del kit IL / Nom du kit IL

In vitro diagnostic medical device De uso diagnóstico in vitro Dispositif médical de diagnostic in vitro	Batch code Identificación número de lote Désignation du lot	Use by Caducidad Utilisable jusqu'à	Temperature limitation Temperatura de Almacenamiento Températures limites de conservation	Consult instructions for use Consultar la metódica Lire le mode d'emploi	Control Control Contrôle	Biological risks Riesgo biológico Risque biologique	Manufacturer Fabricado por Fabricant	Change(s) to section since previous revision La sección ha cambiado desde la revisión anterior La section a changé depuis la révision précédente	Caution, consult accompanying documents Atención, ver instrucciones de uso Attention, voir notice d'instructions
--	---	---	---	--	--------------------------------	---	--	--	--

Printed Insert Sheet: 303418
Revision: R16
Issued: 03/2020
C.O.: 515534

LANGUAGES

ENGLISH
SPANISH
FRENCH

TECHNICAL SPECS

PAPER: White paper, 50-60 g/m² weight, PCS00400073.
SIZE: 11 x 17" (280 x 432 mm.).
PRINT: Front/Back.
PRINT COLOR: Front - Top band Green Pantone 382, all remaining
type in black.
Back - All type in black.
SPECIAL REQUIREMENTS: Folded