



DIAPHARMA I-2581
Catalog # B821108-25

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR HUMAN USE OR DRUG USE.

DIAPHARMA I-2581

COMPOSITION

Each vial contains 25 mg of synthetic thrombin inhibitor and 50 mg mannitol added as a bulking agent.

Chromogenic Substrate –

CHEMISTRY

N^α-dansyl-(p-guanidino)-phenylalanine-piperidide hydrochloride

MW: 559.2

Solubility: > 10 mmol/L

Suitable stock solution: 0.5 mmol/L in H₂O

STORAGE

May be used by the retest date given on the label when store unopened, protected from moisture, in the dark, 2-8°C. Avoid contamination of the reagents by microorganisms.

Diapharma Group, Inc.
8948 Beckett Road
West Chester, OH 45069
Customer Service:
1-800-526-5224
Technical Service:
1-800-447-3846
Fax:
1- 513- 860-9635
E-mail:
info@diapharma.com
Web:
www.diapharma.com