



DIAPHARMA CS – GK

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR HUMAN USE OR DRUG USE.

DIAPHARMA CS – GK

DiaPharma CS – Glandular Kallikrein is a chromogenic peptide substrate for glandular kallikrein

COMPOSITION

Each vial contains chromogenic substrate H-D-Val-Leu-Arg-pNA · 2HCl 25 mg and mannitol 60 mg added as a bulking agent.

Chromogenic Substrate – Glandular Kallikrein

CHEMISTRY

H-D-Val-Leu-Arg-pNA · 2HCl
MW: 579.6

PRINCIPLE

H-D-Val-Leu-Arg-pNA + Enzyme → H-D-Val-Leu-Arg-OH + pNA

STORAGE

May be used by the expiry date given on the label when store unopened, protected from moisture, in the dark, 2-8°C. Avoid contamination of the reagents by microorganisms. Shipment of product does not require cooling during the time of transportation.

DiaPharma Group, Inc.
8948 Beckett Road
West Chester, OH 45069
Customer Service: 1-800-526-5224
Technical Service: 1-800-447-3846
Fax: (513) 860-9635
E-mail: info@diapharma.com
Web: www.diapharma.com

DMR-07-00002 Rev01